

## Perovskite MCPT Run Protocol

<b>Sample:</b>		<b>ID:</b>		<b>Folder:</b>	
<b>Feed:</b>				<b>Date:</b>	

### Pre-run checklist

Gas valves:			Trace heating:		Temperatures:			Cooling temperatures:	
	top	mfc				Actual:	Setpoint:		
N2:			Preheat:		260°C			Peltier	
O2:			GC 1:		270°C			Pump	
alkane:			GC 2:		280°C				
flush:			GC 3:		290°C				
heater:					300°C				
					calibration:				

### Sample & flows:

Sample properties:				Flow mixtures:			
Mass:				Inert			
before					actual	setpoint	
after				O2	10%		
difference				N2	90%		
Height:				Feed			
Height:					%:	actual	setpoint
				alkane			
				+ N2			
				O2			
				N2			

### Flow setpoints [gs/ml] and [ml/min]

	res time	flow	high	low	set high	set low
inert	0.20					
	0.20					
feed	0.25					
	0.30					
calibration:						

### GC

Method:	
calibration:	

### MCPT

Bandwidth		Start	
Points		Stop	
Averaging		Calibration	
Cold Coupling			
	Freq.	log( a+bi )	
TM210			
TM020			